

Strengthen Confidence with High Pressure

BRAVO™

Non Compliant Balloon Dilatation Catheter





| Highest RBP | Low Compliance | Hydrophilic Coating |

SINOMED
Innovation for health

BRAVO™


NC Balloon Dilatation

With brand new balloon material, the BRAVO NC balloon catheter is designed for optimal post-dilatation of coronary stents.



Highest RBP to boost confidence for physicians aiming to optimal dilatation.

on Catheter



Lower compliance to achieve superior stent apposition.



Low crossing profile to deliver superb lesion access.



Selective hydrophilic coating to enhance the crossability without compromising precise placement.

BRAVO™ NC Balloon Dilatation Catheter

Technical Information

Balloon Material	Nylon
Shaft Diameters	2.0F proximal, 2.6F distal
Usable Length	145 cm
Minimum Guiding Catheter	5F (0.056")
Maximum Guide Wire	0.014"

Order Information

Length (mm)	Diameter (mm)									
	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.25	4.50	5.00
6	●	●	●	●	●	●	●	●	●	●
10	●	●	●	●	●	●	●	●	●	●
15	●	●	●	●	●	●	●	●	●	●
20	●	●	●	●	●	●	●	●	●	●

Compliance Table

Pressure (atm)	Diameter (mm)									
	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.25	4.50	5.00
6	2.36	2.49	2.55	2.96	3.09	3.36	3.59	3.85	4.05	4.56
8	2.41	2.58	2.71	3.07	3.24	3.49	3.71	3.98	4.21	4.71
10	2.46	2.69	2.86	3.17	3.37	3.61	3.84	4.11	4.35	4.85
12	2.51	2.75	2.97	3.24	3.51	3.74	4.02	4.25	4.48	4.98
14	2.53	2.78	3.01	3.28	3.55	3.78	4.05	4.29	4.51	5.03
16	2.54	2.80	3.03	3.31	3.59	3.82	4.09	4.32	4.55	5.09
18	2.56	2.82	3.06	3.35	3.64	3.86	4.13	4.36	4.59	5.13
20	2.58	2.84	3.10	3.38	3.69	3.89	4.17	4.40	4.64	5.18
22	2.59	2.85	3.13	3.41	3.73	3.93	4.20	4.43	4.68	5.24
24	2.61	2.87	3.15	3.43	3.76	3.96	4.22			

■ NP ■ RBP

Sino Medical Sciences Technology Inc

2nd Floor, TEDA Biopharm Res, Building B, #5 4th St, TEDA, Tianjin, China

Z : 300457

T : + 86 22 59885298

F : + 86 22 62000060

E: globalbusiness@sinomedical.net

en.sinomedical.net

